



TABLE OF AVAILABLE COIL CONFIGURATIONS AND CORRESPONDING COILS

1H RES 23mm Quad	Mouse Head Volume Coil
1H RES 40mm Quad	Rat Head/Mouse Body Volume Coil
1H RES 60mm Quad	Rat Body Volume Coil
LNA Loop Coil RES 72mm TO	Circular Loop Coils (10/20/30mm)
M.Br. Arr RES 72mm TO	Mouse Brain Array Surface Coil
QuadTransceiver 1H 300W	Old Rat Body Volume Coil from Rapid
R.Br.Arr. RES 72mm	Rat Brain Array Surface Coil
Transceiver 1H 300W	Old Rat Head Volume Coil
Transceiver 1H 300 Coil array	
Transceiver 1H 300W CrossCoil	